



USER MANUAL – EN

IN 30614 Mechanical Folding Sports Wheelchair inSPORTline Brenstark Lightweight



The images in the manual are for illustrative purposes only, the final product finish may slightly vary.

SevenSport s.r.o. reserves the right to make any changes and improvements to its product without prior notice. Visit our website www.insportline.eu where you will find the latest version of the manual.



CONTENTS

SAFETY INSTRUCTIONS.....	3
CAUTIONS AND WARNINGS.....	4
ADJUSTABLE PARTS	4
Notes.....	5
MAINTENANCE	5
SPECIFICATIONS.....	6
ENVIRONMENT PROTECTION	6
TERMS AND CONDITIONS OF WARRANTY, WARRANTY CLAIMS.....	6

SAFETY INSTRUCTIONS

- Read the manual before use and keep it for future reference.
- Failure to follow the instructions and warnings may result in injury, property damage or death.
- The stability of the wheelchair is affected by the condition of the front and rear wheels, as well as the user.
- Do not drive uphill or downhill with an incline greater than 6°.
- Do not drive on public motor roads.
- Do not drive over obstacles such as curbs or gaps. The wheelchair may tip over.
- Do not climb into the wheelchair.
- Do not use the footrests as a platform. When dismounting, make sure the footrests are folded.
- Always use the wheel handles when pushing yourself. Do not use the tires.
- Do not use other unapproved accessories or attachments.
- The temperature of the wheelchair may increase if the wheelchair is exposed to external heat sources (sunlight).
- Use the wheelchair only for the purpose for which it was manufactured.
- The seller or manufacturer is not responsible for injuries or property damage due to improper use of the wheelchair.
- The user and assistant must be instructed in proper use.
- Keep limbs and body out of the wheelchair.
- Keep the center of gravity in the center of the seat. Always sit upright and do not lean forward.
- Regularly check the tightness of all joints and adjustable parts. Never lift the wheelchair by the adjustable parts.
- Regularly check the condition of the wheelchair, do not use a worn or damaged wheelchair.
- When driving over obstacles or edges, be careful not to tip the wheelchair over. Never let the wheelchair fall under its weight.
- When using a lift, pay attention to the load capacity and type of lift.
- Do not lift the wheelchair if a person is sitting in it. There is a risk of injury.
- If it is necessary to transfer a person from the wheelchair to another surface, try to reduce the distance between the two places as much as possible.
- When transferring a person from the wheelchair, remember to properly secure the wheelchair.
- Keep limbs and body parts away from moving parts.
- The wheelchair is designed to carry one person at a time.
- Only carry out any repairs, maintenance or cleaning when the wheelchair is empty.
- Do not overtighten adjustable parts and joints. This can damage the frames and threads.
- Adjustable parts must be able to be adjusted without the need for tools.
- Only use identical original spare parts for repairs.
- If sharp edges appear on the wheelchair, stop using the wheelchair immediately until the damage has been repaired.
- Load capacity: 100 kg

CAUTIONS AND WARNINGS

		
Do not lean over the edge of the cart unless necessary. If necessary, lock the wheels.	Do not lean back, you may change your centre of gravity and tip over. se.	When driving over an obstacle, tilt the wheelchair by stepping on the frame.

ADJUSTABLE PARTS

Press the end of the rear wheel axle to remove or install the rear wheels.	
Use the Allen screw to adjust the height of the rear safety wheels. Tighten the bolts after adjustment.	
Height-adjustable footrest. Adjust the height using the Allen screws. Then tighten the screws.	

Handbrakes on both sides.



Notes

Check the position of the safety wheels before use. Adjust them so that they do not hinder movement before pushing.



Always use a seat belt.



Check the security of the adjustable parts.



MAINTENANCE

- For cleaning, use only mild detergents and lukewarm water. Then dry thoroughly. Do not use aggressive detergents, bleaches, etc.
- Keep the wheels clean. Dirt can slow down or completely stop their movement.
- Regularly lubricate all moving parts.
- Sand, salt water or salt can damage metal parts, bearings and other moving parts, and may cause corrosion of metal parts.
- Regularly check the condition of all moving joints and contact the service department if necessary.

SPECIFICATIONS

Wheelchair part	Dimension
Total length	85 cm
Total width	58 cm
Total height	80 cm
Unfolded width	35 cm
Folded width	36 cm
Seat depth	36 cm
Seat to ground height	48 cm
Footrest to ground height	5-10 cm
Backrest height	27 cm
Front wheel diameter	10 cm
Rear wheel diameter	60 cm
Max. load capacity	100 kg
Net weight	10 kg

ENVIRONMENT PROTECTION

After the product lifespan expired or if the possible repairing is uneconomic, dispose it according to the local laws and environmentally friendly in the nearest scrapyard.

By proper disposal you will protect the environment and natural sources. Moreover, you can help protect human health. If you are not sure in correct disposing, ask local authorities to avoid law violation or sanctions.

Don't put the batteries among house waste but hand them in to the recycling place.

TERMS AND CONDITIONS OF WARRANTY, WARRANTY CLAIMS

General Conditions of Warranty and Definition of Terms

All Warranty Conditions stated here under determine Warranty Coverage and Warranty Claim Procedure. Conditions of Warranty and Warranty Claims are governed by Act No. 89/2012 Coll. Civil Code, and Act No. 634/1992 Coll., Consumer Protection, as amended, also in cases that are not specified by these Warranty rules.

The seller is SEVEN SPORT s.r.o. with its registered office in Kutnohorská 531, Štěrboholy Retail Park, 109 00, Praha 10, ČR, Company Registration Number: 26847264, registered in the Trade Register at Regional Court in Prague, Section C, Insert No. 116888.

According to valid legal regulations it depends whether the Buyer is the End Customer or not.

"The Buyer who is the End Customer" or simply the "End Customer" is the legal entity that does not conclude and execute the Contract in order to run or promote his own trade or business activities.

"The Buyer who is not the End Customer" is a Businessman that buys Goods or uses services for the purpose of using the Goods or services for his own business activities. The Buyer conforms to the General Purchase Agreement and business conditions.

These Conditions of Warranty and Warranty Claims are an integral part of every Purchase Agreement made between the Seller and the Buyer. All Warranty Conditions are valid and binding, unless otherwise specified in the Purchase Agreement, in the Amendment to this Contract or in another written agreement.

Warranty Conditions

Warranty Period

The Seller provides the Buyer a 24 months Warranty for Goods Quality, unless otherwise specified in the Certificate of Warranty, Invoice, Bill of Delivery or other documents related to the Goods. The legal warranty period provided to the Consumer is not affected.

By the Warranty for Goods Quality, the Seller guarantees that the delivered Goods shall be, for a certain period of time, suitable for regular or contracted use, and that the Goods shall maintain its regular or contracted features.

Batteries

6-month battery warranty – we guarantee that battery's nominal capacity does not fall below 70% of its total capacity within 6 months of the product's sale.

The Warranty does not cover defects resulting from (if applicable):

- User's fault, i.e. product damage caused by unqualified repair work, improper assembly, insufficient insertion of seat post into frame, insufficient tightening of pedals and cranks
- Improper maintenance
- Mechanical damages
- Regular use (e.g. wearing out of rubber and plastic parts, moving mechanisms, joints, wear of brake pads/blocks, chain, tires, cassette/multi wheel etc.)
- Unavoidable event, natural disaster
- Adjustments made by unqualified person
- Improper maintenance, improper placement, damages caused by low or high temperature, water, inappropriate pressure, shocks, intentional changes in design or construction etc.

Warranty Claim Procedure

The Buyer is obliged to check the Goods delivered by the Seller immediately after taking the responsibility for the Goods and its damages, i.e. immediately after its delivery. The Buyer must check the Goods so that he discovers all the defects that can be discovered by such check.

When making a Warranty Claim the Buyer is obliged, on request of the Seller, to prove the purchase and validity of the claim by the Invoice or Bill of Delivery that includes the product's serial number, or eventually by the documents without the serial number. If the Buyer does not prove the validity of the Warranty Claim by these documents, the Seller has the right to reject the Warranty Claim.

If the Buyer gives notice of a defect that is not covered by the Warranty (e.g. in the case that the Warranty Conditions were not fulfilled or in the case of reporting the defect by mistake etc.), the Seller is eligible to require a compensation for all the costs arising from the repair. The cost shall be calculated according to the valid price list of services and transport costs.

If the Seller finds out (by testing) that the product is not damaged, the Warranty Claim is not accepted. The Seller reserves the right to claim a compensation for costs arising from the false Warranty Claim.

In case the Buyer makes a claim about the Goods that is legally covered by the Warranty provided by the Seller, the Seller shall fix the reported defects by means of repair or by the exchange of the damaged part or product for a new one. Based on the agreement of the Buyer, the Seller has the right to exchange the defected Goods for a fully compatible Goods of the same or better technical characteristics. The Seller is entitled to choose the form of the Warranty Claim Procedures described in this paragraph.

The Seller shall settle the Warranty Claim within 30 days after the delivery of the defective Goods, unless a longer period has been agreed upon. The day when the repaired or exchanged Goods is handed over to the Buyer is considered to be the day of the Warranty Claim settlement. When the Seller is not able to settle the Warranty Claim within the agreed period due to the specific nature of the Goods defect, he and the Buyer shall make an agreement about an alternative solution. In case such agreement is not made, the Seller is obliged to provide the Buyer with a financial compensation in the form of a refund.

CZ
SEVEN SPORT s.r.o.

Registered Office: Kutnohorská 531, Štěrboholy Retail Park,
109 00, Praha 10, ČR

Headquaters: Dělnická 957, Vítkov, 749 01
Warranty & Service: Čermenská 486, Vítkov 749 01

CRN: 26847264
VAT ID: CZ26847264

Phone: +420 556 300 970
E-mail: eshop@insportline.cz
reklamace@insportline.cz
servis@insportline.cz

Web: www.inSPORTline.cz

SK
Stores inSPORTline SK, s.r.o.

Headquaters, warranty & service center: Električná 6471,
Trenčín 911 01, SK

CRN: 46259317
VAT ID: SK2023299729

Phone: +421(0) 917 700 098
E-mail: objednavky@insportline.sk
reklamacie@insportline.sk
servis@insportline.sk

Web: www.inSPORTline.sk

About shipping

