



**USER MANUAL – EN**

**IN 18010 Feet massager inSPORTline Otterchill**



**CONTENTS**

SAFETY INSTRUCTIONS..... 3  
USE ..... 3  
BATTERY REPLACEMENT ..... 4  
ALTERNATIVE USE..... 4  
TERMS AND CONDITIONS OF WARRANTY, WARRANTY CLAIMS..... 5

## **SAFETY INSTRUCTIONS**

- Read the manual carefully.
- Always unplug the product after use.
- If the device come in touch with the water, do not touch it and unplug it immediately.
- Do not use while taking a shower or bathing.
- Do not store in humid places, in or near water.
- Do not drop in liquid.
- Do not leave unattended when plugged in. Always unplug the product after use.
- Do not use under blanket or pillow. Excessive heating can cause fire, electric shock or injuries.
- Supervision is necessary when using the product near children or people with reduce mental or physical capabilities.
- Use only for intended purposes. Do not modify the product.
- Never try to use the product when the power cord is damaged or wet. If the product is damaged ask your supplier or expert for repair.
- Do not pull the power cord.
- Keep away the power cord from sources of heat.
- Always keep the ventilation free.
- Never drop any items into the ventilation.
- Do not use outside.
- Turn off all functions and then unplug from socket.
- Do not use the heating function when sensitive to heat. Keep away from kids and persons with reduce mental or physical capabilities.
- Use the massager only when sitting.
- Do not run the massager for longer than 1 hour. It is recommended to use the massager in intervals of 20 minutes.
- People suffering from: chronic or severe skin disease, hearth problems, physical injury, fever, stroke or pregnant women and mentally ill people should not use the product.
- Do not use under the influence of drugs and alcohol.
- Do not use within one hour before and after meal.
- To clean the product use damp cloth. Do not use alcohol or water to clean the product.
- Do not repair or modify the product on your own.
- Do not stand or sit on the massager.
- To eliminate skin irritation wear socks during the massage.
- Do not drop the massager.
- Do not inset foreign objects into the air vent.


## **USE**

- Plug to socket.
- Display will show OF which means that the massager is ready and in standby mode. Start the massage by pressing ON/OFF button. Display will show OO. Massager will automatically at 1 vibration level.

- Adjust the speed by pressing + and – button. There are 15 vibration levels. The display will switch to calories after the speed is selected.
- To turn on the heat function – press the FIR button on the device and Infrared on the remote controller. Heat function can be turned on independently.
- Timer can be set to 10, 20, 30 minutes. Massage will stop automatically after selected time.
- For auto pre-set massage press the SELECT button. Press + and – to go back to manual mode.
- 4 pre-sets modes:
  - Vitality mode
  - Relaxation mode
  - Reflexology mode
  - Slimming mode
- When using pre-set massage, you can select time.
- Calories count will be displayed after few seconds.
- The massage will start he low level, after that you can adjust the level to your requirements. Do not use for more than 30minutes.
- Remove the plug form socket after use.




## BATTERY REPLACEMENT

- Battery is located on the back of the wireless controller.
- Press the small clip and pull out the drawer battery compartment.
- Use “LITHIUM” battery CR2025.



## ALTERNATIVE USE

The massager can be used on other parts of body. The vibration helps to relax muscles and reduce stiffens when used on calves, thighs, back or waist. It can be used regularly on abdominal and thigh areas.

	<p>Place your feet on the center of the pressure pad.</p>
	<p>Lean your upper back and shoulder on the pressure pad to reduce pain on neck, arms or shoulders.</p>
	<p>Sit comfortably and place the center of the massager's pressure pad under your thighs. (Wear long trousers.)</p>

## TERMS AND CONDITIONS OF WARRANTY, WARRANTY CLAIMS

### General Conditions of Warranty and Definition of Terms

All Warranty Conditions stated hereunder determine Warranty Coverage and Warranty Claim Procedure. Conditions of Warranty and Warranty Claims are governed by Act No. 40/1964 Coll. Civil Code, Act No. 513/1991 Coll., Commercial Code, and Act No. 634/1992 Coll., Consumer Protection Act, as amended, also in cases that are not specified by these Warranty rules.

The seller is SEVEN SPORT s.r.o. with its registered office in Borivojova Street 35/878, Prague 13000, Company Registration Number: 26847264, registered in the Trade Register at Regional Court in Prague, Section C, Insert No. 116888.

According to valid legal regulations it depends whether the Buyer is the End Customer or not.

“The Buyer who is the End Customer” or simply the “End Customer” is the legal entity that does not conclude and execute the Contract in order to run or promote his own trade or business activities.

“The Buyer who is not the End Customer” is a Businessman that buys Goods or uses services for the purpose of using the Goods or services for his own business activities. The Buyer conforms to the General Purchase Agreement and business conditions to the extent specified in the Commercial Code.

These Conditions of Warranty and Warranty Claims are an integral part of every Purchase Agreement made between the Seller and the Buyer. All Warranty Conditions are valid and binding, unless

otherwise specified in the Purchase Agreement, in the Amendment to this Contract or in another written agreement.

### **Warranty Conditions**

#### **Warranty Period**

The Seller provides the Buyer a 24 months Warranty for Goods Quality, unless otherwise specified in the Certificate of Warranty, Invoice, Bill of Delivery or other documents related to the Goods. The legal warranty period provided to the Consumer is not affected.

By the Warranty for Goods Quality, the Seller guarantees that the delivered Goods shall be, for a certain period of time, suitable for regular or contracted use, and that the Goods shall maintain its regular or contracted features.

#### **The Warranty does not cover defects resulting from (if applicable):**

- User's fault, i.e. product damage caused by unqualified repair work, improper assembly, insufficient insertion of seat post into frame, insufficient tightening of pedals and cranks
- Improper maintenance
- Mechanical damages
- Regular use (e.g. wearing out of rubber and plastic parts, moving mechanisms, joints etc.)
- Unavoidable event, natural disaster
- Adjustments made by unqualified person
- Improper maintenance, improper placement, damages caused by low or high temperature, water, inappropriate pressure, shocks, intentional changes in design or construction etc.

### **Warranty Claim Procedure**

The Buyer is obliged to check the Goods delivered by the Seller immediately after taking the responsibility for the Goods and its damages, i.e. immediately after its delivery. The Buyer must check the Goods so that he discovers all the defects that can be discovered by such check.

When making a Warranty Claim the Buyer is obliged, on request of the Seller, to prove the purchase and validity of the claim by the Invoice or Bill of Delivery that includes the product's serial number, or eventually by the documents without the serial number. If the Buyer does not prove the validity of the Warranty Claim by these documents, the Seller has the right to reject the Warranty Claim.

If the Buyer gives notice of a defect that is not covered by the Warranty (e.g. in the case that the Warranty Conditions were not fulfilled or in the case of reporting the defect by mistake etc.), the Seller is eligible to require a compensation for all the costs arising from the repair. The cost shall be calculated according to the valid price list of services and transport costs.

If the Seller finds out (by testing) that the product is not damaged, the Warranty Claim is not accepted. The Seller reserves the right to claim a compensation for costs arising from the false Warranty Claim.

In case the Buyer makes a claim about the Goods that is legally covered by the Warranty provided by the Seller, the Seller shall fix the reported defects by means of repair or by the exchange of the damaged part or product for a new one. Based on the agreement of the Buyer, the Seller has the right to exchange the defected Goods for a fully compatible Goods of the same or better technical characteristics. The Seller is entitled to choose the form of the Warranty Claim Procedures described in this paragraph.

The Seller shall settle the Warranty Claim within 30 days after the delivery of the defective Goods, unless a longer period has been agreed upon. The day when the repaired or exchanged Goods is handed over to the Buyer is considered to be the day of the Warranty Claim settlement. When the Seller is not able to settle the Warranty Claim within the agreed period due to the specific nature of the Goods defect, he and the Buyer shall make an agreement about an alternative solution. In case such agreement is not made, the Seller is obliged to provide the Buyer with a financial compensation in the form of a refund.



**SEVEN SPORT s.r.o.**

Registered Office: Borivojova 35/878, 130 00 Praha 3, Czech Republic  
Headquarters: Delnicka 957, 749 01 Vitkov, Czech Republic  
Warranty & Service Centre: Cermenska 486, 749 01 Vitkov, Czech Republic  
CRN: 26847264  
VAT ID: CZ26847264  
Phone: +420 556 300 970  
E-mail: eshop@insportline.cz  
reklamace@insportline.cz  
servis@insportline.cz  
Web: www.insportline.cz

---

**SK**

**INSPORTLINE s.r.o.**

Headquarters, Warranty & Service centre: Elektricka 6471, 911 01 Trencin, Slovakia  
CRN: 36311723  
VAT ID: SK2020177082  
Phone: +421(0)326 526 701  
E-mail: objednavky@insportline.sk  
reklamacie@insportline.sk  
servis@insportline.sk  
Web: www.insportline.sk

Date of Sale:

Stamp and Signature of Seller: