



USER MANUAL – EN

IN 1992 Massage collar inSPORTline S002



CONTENTS

SAFETY INSTRUCTIONS..... 3
USE 4
PRODUCT DESCRIPTION 5
CONTROLS..... 6
TROUBLESHOOTING 6
TECHNICAL PARAMETERS 6
ENVIRONMENT PROTECTION 7
TERMS AND CONDITIONS OF WARRANTY, WARRANTY CLAIMS..... 7

SAFETY INSTRUCTIONS


- Read the manual carefully and keep it for future reference.
- Not suitable for children or pregnant women.
- Do not use while driving.
- Do not use if you suffer from the following health problems: heat sensitivity, circulation problems, varicose veins, bruising, skin irritation, inflammation of muscles or tissues, inflammation of the veins or open wounds
- If you experience pain or discomfort, turn off the device immediately and contact a physician.
- Do not use the product if you are using a pacemaker or other similar medical implants.
- Do not use if you have used cream or oil on your skin.
- Turn off the product and unplug it if not in use.
- Do not sleep during use.
- Do not transport the product or pull by the power cord.
- Clean the product if it is not plugged in. Use only a damp cloth for cleaning. Then dry the product immediately. Do not connect until the device is dry.
- Do not cover any blanket or lie on a pillow. Accumulated heat can damage the product or injure the user.
- Do not use in places with high humidity or near water. Do not use the product if it is not completely dry.
- Do not expose to any liquids.
- The product is for home use only.
- Keep out of reach of children.
- Do not leave the product unattended when it is turned on.
- Do not pull or leave the power cord in places where there is a risk of tripping.
- Do not leave near a heat source.
- Do not use the product if it is damaged.
- Check the product and power cord regularly for damage. If it needs a repair, contact your dealer, authorized service or specialist.
- Use only original spare parts for repair.
- If the product is damaged, stop it immediately.
- Never use the product if wet hands.
- Never open or repair the product at home.
- Do not use at bedtime, massage may interfere with sleep.
- Keep all packing materials out of the reach of children.

USE

Adjust the massage intensity to your needs.




Neck&Shoulder




Neck



Back



Waist



Abdomen



Thighs



Legs






Foot







Calf

PRODUCT DESCRIPTION

<p>A</p>  <p>Diagram of the main unit showing callouts: 1 - Leather cover, 2 - Massage head, 3 - Power switch, 4 - Handle, 5 - DC connector.</p>	<p>A – Main unit</p> <p>1 – Leather cover</p> <p>2 – Massage head</p> <p>3 – Power switch</p> <p>4 – Handle</p> <p>5 – DC connector</p>
<p>B</p>  <p>Diagram of the adapter showing callouts: 6 - Plug, 7 - Cord, 8 - Appliance plug.</p>	<p>B – Adapter</p> <p>6 – Plug</p> <p>7 – Cord</p> <p>8 – Appliance plug</p>
<p>C</p>  <p>Diagram of the car cigarette lighter cord.</p>	<p>C – Car cigarette lighter</p>

CONTROLS

	<p>Turn on / off Turn on / off the device. The device will automatically turn off after 15 minutes.</p>
	<p>Speed Adjust the speed of massage</p>
	<p>Direction Adjust the direction of massage.</p>
	<p>Heating Turn on / off</p>

TROUBLESHOOTING

PROBLEM	SOLUTION
Unable to start massage	Make sure that the plug is in the socket. Confirm that the switch is on.
Suddenly stops working during use	<p>Massage will automatically turn off after 15 minutes.</p> <p>Make sure that the plug is in the socket.</p> <p>Check the plug and power cord. If they are damaged, stop using it immediately.</p> <p>Product has overheat protection. Turn it off for at least one hour.</p> <p>If the problems persist, contact your local dealer.</p>

TECHNICAL PARAMETERS

Voltage: Dc 12 V

Power: 24 W

Frequency: 50 Hz

Product features or design might be changed without notice.

ENVIRONMENT PROTECTION

After the product lifespan expired or if the possible repairing is uneconomic, dispose it according to the local laws and environmentally friendly in the nearest scrapyards.

By proper disposal you will protect the environment and natural sources. Moreover, you can help protect human health. If you are not sure in correct disposing, ask local authorities to avoid law violation or sanctions

TERMS AND CONDITIONS OF WARRANTY, WARRANTY CLAIMS

General Conditions of Warranty and Definition of Terms

All Warranty Conditions stated hereunder determine Warranty Coverage and Warranty Claim Procedure. Conditions of Warranty and Warranty Claims are governed by Act No. 89/2012 Coll. Civil Code, and Act No. 634/1992 Coll., Consumer Protection, as amended, also in cases that are not specified by these Warranty rules.

The seller is SEVEN SPORT s.r.o. with its registered office in Borivojova Street 35/878, Prague 13000, Company Registration Number: 26847264, registered in the Trade Register at Regional Court in Prague, Section C, Insert No. 116888.

According to valid legal regulations it depends whether the Buyer is the End Customer or not.

“The Buyer who is the End Customer” or simply the “End Customer” is the legal entity that does not conclude and execute the Contract in order to run or promote his own trade or business activities.

“The Buyer who is not the End Customer” is a Businessman that buys Goods or uses services for the purpose of using the Goods or services for his own business activities. The Buyer conforms to the General Purchase Agreement and business conditions.

These Conditions of Warranty and Warranty Claims are an integral part of every Purchase Agreement made between the Seller and the Buyer. All Warranty Conditions are valid and binding, unless otherwise specified in the Purchase Agreement, in the Amendment to this Contract or in another written agreement.

Warranty Conditions

Warranty Period

The Seller provides the Buyer a 24 months Warranty for Goods Quality, unless otherwise specified in the Certificate of Warranty, Invoice, Bill of Delivery or other documents related to the Goods. The legal warranty period provided to the Consumer is not affected.

By the Warranty for Goods Quality, the Seller guarantees that the delivered Goods shall be, for a certain period of time, suitable for regular or contracted use, and that the Goods shall maintain its regular or contracted features.

The Warranty does not cover defects resulting from (if applicable):

- User's fault, i.e. product damage caused by unqualified repair work, improper assembly, insufficient insertion of seat post into frame, insufficient tightening of pedals and cranks
- Improper maintenance
- Mechanical damages
- Regular use (e.g. wearing out of rubber and plastic parts, moving mechanisms, joints, wear of brake pads/blocks, chain, tires, cassette/multi wheel etc.)
- Unavoidable event, natural disaster
- Adjustments made by unqualified person

- Improper maintenance, improper placement, damages caused by low or high temperature, water, inappropriate pressure, shocks, intentional changes in design or construction etc.

Warranty Claim Procedure

The Buyer is obliged to check the Goods delivered by the Seller immediately after taking the responsibility for the Goods and its damages, i.e. immediately after its delivery. The Buyer must check the Goods so that he discovers all the defects that can be discovered by such check.

When making a Warranty Claim the Buyer is obliged, on request of the Seller, to prove the purchase and validity of the claim by the Invoice or Bill of Delivery that includes the product's serial number, or eventually by the documents without the serial number. If the Buyer does not prove the validity of the Warranty Claim by these documents, the Seller has the right to reject the Warranty Claim.

If the Buyer gives notice of a defect that is not covered by the Warranty (e.g. in the case that the Warranty Conditions were not fulfilled or in the case of reporting the defect by mistake etc.), the Seller is eligible to require a compensation for all the costs arising from the repair. The cost shall be calculated according to the valid price list of services and transport costs.

If the Seller finds out (by testing) that the product is not damaged, the Warranty Claim is not accepted. The Seller reserves the right to claim a compensation for costs arising from the false Warranty Claim.

In case the Buyer makes a claim about the Goods that is legally covered by the Warranty provided by the Seller, the Seller shall fix the reported defects by means of repair or by the exchange of the damaged part or product for a new one. Based on the agreement of the Buyer, the Seller has the right to exchange the defected Goods for a fully compatible Goods of the same or better technical characteristics. The Seller is entitled to choose the form of the Warranty Claim Procedures described in this paragraph.

The Seller shall settle the Warranty Claim within 30 days after the delivery of the defective Goods, unless a longer period has been agreed upon. The day when the repaired or exchanged Goods is handed over to the Buyer is considered to be the day of the Warranty Claim settlement. When the Seller is not able to settle the Warranty Claim within the agreed period due to the specific nature of the Goods defect, he and the Buyer shall make an agreement about an alternative solution. In case such agreement is not made, the Seller is obliged to provide the Buyer with a financial compensation in the form of a refund.

CZ

SEVEN SPORT s.r.o.

Registered Office: Strakonická 1151/2c, Praha 5, 150 00, ČR
 Headquarters: Dělnická 957, Vítkov, 749 01
 Warranty & Service: Čermenská 486, Vítkov 749 01

CRN: 26847264
 VAT ID: CZ26847264

Phone: +420 556 300 970
 E-mail: eshop@insportline.cz
 reklamace@insportline.cz
 servis@insportline.cz

Web: www.inSPORTline.cz

SK

inSPORTline s.r.o.

Headquarters, warranty & service center: Električná 6471, Trenčín 911 01, SK

CRN: 36311723
 VAT ID: SK2020177082

Phone: +421(0)326 526 701
 E-mail: objednavky@insportline.cz
 reklamacie@insportline.cz
 servis@insportline.cz

Web: www.inSPORTline.sk