

USER MANUAL – EN

IN 23199 Ultrasonic massage device inSPORTline Valeri


CONTENTS

SAFETY INSTRUCTIONS	. 3
PRODUCT DESCRIPTION	. 3
USE - ULTRASONIC PROCEDURE	. 4
TIPS FOR PLACING ELECTRODES	. 6
AFTER USE	. 7
ENVIRONMENT PROTECTION	. 7
TERMS AND CONDITIONS OF WARRANTY, WARRANTY CLAIMS	. 7

SAFETY INSTRUCTIONS


- Read the manual before use and keep it for future reference.
- Do not use the device if the power cord is hot or the plug is damaged.
- Keep away from water and moisture.
- The device is not a toy, keep out of reach of children and pets.
- Do not move the electrodes or allow another person to touch the device during use.
- Do not place the electrodes on the eyelids, forehead, chest or near the heart.
- Not suitable for: pregnant women, people suffering from infection, dermatological diseases, people using pacemakers or other electrical devices, etc.
- Always set the device to the lowest intensity before starting or switching off.
- WARNING: The device is not a substitute for medical or therapeutic care.

PRODUCT DESCRIPTION



USE - ULTRASONIC PROCEDURE

- 1. Connect the device to the power supply.
- 2. Start the device, the red indicator lights up.




- 3. Press the SONIC button. The device starts at the lowest intensity. To increase the intensity, press the SONIC button repeatedly. For more sensitive areas, we recommend low intensity.
- 4. For used areas, apply the conductive gel or cream so that the skin does not dry out.



- 5. Apply cream on the head.
- 6. Massage the skin on the desired areas with circular movements, see the pictures below



WARNING: Do not use the device on the neck, eyelids, eyes, joints or bones.



USE – EMS

- 1. Clean the necessary part of the body from grease, cream, etc. with a damp towel.
- 2. Connect the electrodes to the cables and then place the positive and negative electrodes on different parts of the body.







- 3. Set the intensity to the lowest value.
- 4. Connect the cables to the device.

- 5. Start the device (red LED lights up).
- 6. Choose the type of massage (Tapping, massage, knead, scrapping, slimming)
- 7. Press the EMS button.
- 8. Set the intensity to one of 5 values.
- 9. The device switches off automatically after 10 minutes of operation.
- 10. We recommend starting at 10 minutes and possibly increasing to a maximum of 30 minutes.
- 11. When done, turn off the selected mode and set the intensity to the lowest value.




TIPS FOR PLACING ELECTRODES

WARNING: The following electrode positions are for reference only and are not for medical or therapeutical purposes.

Acupuncture points and types of massages are described below:

Tapping




Scrapping



Skin beauty: (44): Feiyu (45): Gangyu (5): Roumen

Slimming



Weight reduction:
(42): Dazhui
(52): Yangguan
(8): Zhongyuan
(18): Guanyuan
(22): Juliao
(47): Shenyu
(62): Chengfu

HEATING

To start heating, start the device and press the INFRARED button. After a minute, the device warms up. Do not look at the head when heating.

AFTER USE

- 1. Make sure all modes are turned off.
- 2. Turn off the device.
- 3. Disconnect the electrodes and store them safely.
- 4. Wipe the massage head and electrodes with a damp cloth.
- 5. Do not use solvents or aggressive cleaning agents for cleaning.
- 6. Store the device in the box.

ENVIRONMENT PROTECTION

After the product lifespan expired or if the possible repairing is uneconomic, dispose it according to the local laws and environmentally friendly in the nearest scrapyard.

By proper disposal you will protect the environment and natural sources. Moreover, you can help protect human health. If you are not sure in correct disposing, ask local authorities to avoid law violation or sanctions.

Don't put the batteries among house waste but hand them in to the recycling place.

TERMS AND CONDITIONS OF WARRANTY, WARRANTY CLAIMS

General Conditions of Warranty and Definition of Terms

All Warranty Conditions stated here under determine Warranty Coverage and Warranty Claim Procedure. Conditions of Warranty and Warranty Claims are governed by Act No. 89/2012 Coll. Civil Code, and Act No. 634/1992 Coll., Consumer Protection, as amended, also in cases that are not specified by these Warranty rules.

The seller is SEVEN SPORT s.r.o. with its registered office in Strakonická street 1151/2c, Prague 150 00, Company Registration Number: 26847264, registered in the Trade Register at Regional Court in Prague, Section C, Insert No. 116888.

According to valid legal regulations it depends whether the Buyer is the End Customer or not.

"The Buyer who is the End Customer" or simply the "End Customer" is the legal entity that does not conclude and execute the Contract in order to run or promote his own trade or business activities.

"The Buyer who is not the End Customer" is a Businessman that buys Goods or uses services for the purpose of using the Goods or services for his own business activities. The Buyer conforms to the General Purchase Agreement and business conditions.

These Conditions of Warranty and Warranty Claims are an integral part of every Purchase Agreement made between the Seller and the Buyer. All Warranty Conditions are valid and binding, unless otherwise specified in the Purchase Agreement, in the Amendment to this Contract or in another written agreement.

Warranty Conditions

Warranty Period

The Seller provides the Buyer a 24 months Warranty for Goods Quality, unless otherwise specified in the Certificate of Warranty, Invoice, Bill of Delivery or other documents related to the Goods. The legal warranty period provided to the Consumer is not affected.

By the Warranty for Goods Quality, the Seller guarantees that the delivered Goods shall be, for a certain period of time, suitable for regular or contracted use, and that the Goods shall maintain its regular or contracted features.

Batteries

6-month battery warranty – we guarantee that battery's nominal capacity does not fall below 70% of its total capacity within 6 months of the product's sale.

The Warranty does not cover defects resulting from (if applicable):

- User's fault, i.e. product damage caused by unqualified repair work, improper assembly, insufficient insertion of seat post into frame, insufficient tightening of pedals and cranks
- Improper maintenance
- Mechanical damages
- Regular use (e.g. wearing out of rubber and plastic parts, moving mechanisms, joints, wear of brake pads/blocks, chain, tires, cassette/multi wheel etc.)
- Unavoidable event, natural disaster
- Adjustments made by unqualified person
- Improper maintenance, improper placement, damages caused by low or high temperature, water, inappropriate pressure, shocks, intentional changes in design or construction etc.

Warranty Claim Procedure

The Buyer is obliged to check the Goods delivered by the Seller immediately after taking the responsibility for the Goods and its damages, i.e. immediately after its delivery. The Buyer must check the Goods so that he discovers all the defects that can be discovered by such check.

When making a Warranty Claim the Buyer is obliged, on request of the Seller, to prove the purchase and validity of the claim by the Invoice or Bill of Delivery that includes the product's serial number, or eventually by the documents without the serial number. If the Buyer does not prove the validity of the Warranty Claim by these documents, the Seller has the right to reject the Warranty Claim. If the Buyer gives notice of a defect that is not covered by the Warranty (e.g. in the case that the Warranty Conditions were not fulfilled or in the case of reporting the defect by mistake etc.), the Seller is eligible to require a compensation for all the costs arising from the repair. The cost shall be calculated according to the valid price list of services and transport costs.

If the Seller finds out (by testing) that the product is not damaged, the Warranty Claim is not accepted. The Seller reserves the right to claim a compensation for costs arising from the false Warranty Claim.

In case the Buyer makes a claim about the Goods that is legally covered by the Warranty provided by the Seller, the Seller shall fix the reported defects by means of repair or by the exchange of the damaged part or product for a new one. Based on the agreement of the Buyer, the Seller has the right to exchange the defected Goods for a fully compatible Goods of the same or better technical characteristics. The Seller is entitled to choose the form of the Warranty Claim Procedures described in this paragraph.

The Seller shall settle the Warranty Claim within 30 days after the delivery of the defective Goods, unless a longer period has been agreed upon. The day when the repaired or exchanged Goods is handed over to the Buyer is considered to be the day of the Warranty Claim settlement. When the Seller is not able to settle the Warranty Claim within the agreed period due to the specific nature of the Goods defect, he and the Buyer shall make an agreement about an alternative solution. In case such agreement is not made, the Seller is obliged to provide the Buyer with a financial compensation in the form of a refund.

CZ SEVEN SPORT s.r.o.

Registered Office: ČR Headquaters: Warranty & Service:	Strakonická 1151/2c, Praha 5, 150 00,
	Dělnická 957, Vítkov, 749 01 Čermenská 486, Vítkov 749 01
CRN: VAT ID:	26847264 CZ26847264
Phone: E-mail:	+420 556 300 970 eshop@insportline.cz reklamace@insportline.cz servis@insportline.cz
Web:	www.inSPORTline.cz
SK inSPORTline s.r.o.	

Headquaters, warranty & service center: Električná 6471, Trenčín 911 01, SK

CRN: VAT ID:	36311723 SK2020177082
Phone: E-mail:	+421(0)326 526 701 objednavky@insportline.sk reklamacie@insportline.sk servis@insportline.sk
Web:	www.inSPORTline.sk

About shipping

